

Commuter Companion



Travel Pouch



Use the **Sulcabrush** while reading, watching TV or travelling.

The **Sulcabrush** and the **POCKET Sulcabrush** contribute significantly to improving patients' home care techniques.

For patients who have difficulty flossing!



The **Sulcabrush** and the **POCKET Sulcabrush** are very effective for cleaning around orthodontic bands, crowns, bridges and implants.

**FLIP THE TIP**  
on the **POCKET Sulcabrush**  
to access all areas of the mouth.



Remove and replace tips when they look like this: >



**Sulcabrush** Inc. 21 Bradwick Drive, #12  
Concord, Ontario, Canada  
L4K 1K6

Printed in Canada

Sulcabrush® - U.S. Pat. #4,679,272 / Can. Pat. # 1,247,311

# Sulcabrush®

For patients who have difficulty flossing!

For Healthy Gums!



Travel Pouch



[www.sulcabrush.com/pouch](http://www.sulcabrush.com/pouch)



1-800-387-8777

## Instructions



- < 1. Place the end marked **outside** on the edge of the gumline (photo 1). Continuously follow the wavy gumline, adding a circular motion between the teeth.



- < 2. Continue to apply firm pressure while following the wavy gumline on the **inside** of the mouth, using the end marked inside (photo 2).



- < 3. The **inside** angle is also very effective for reaching the gum behind the last tooth (photo 3).

The **Sulcabrush**<sup>®</sup> was proven better than dental floss in some parts and equal to dental floss in other parts of the teeth and gums for the removal of plaque and the reduction of gingivitis\* (gum infection).

### Studies performed by:

**Sulcabrush**<sup>®</sup> Product Line

Columbia University\*  
School of Dental  
and Oral Surgery

University at Buffalo\*  
School of Dental Medicine

Indiana University\*  
School of Dentistry



[www.sulcabrush.com/pouch](http://www.sulcabrush.com/pouch)

# POCKET™ **Sulcabrush**<sup>®</sup>

Get a grip on healthy Pockets™

For patients  
who have difficulty  
flossing!



If the **Sulcabrush**<sup>®</sup> and the **POCKET™ Sulcabrush**<sup>®</sup> are used daily as directed, they should help stop the gums from bleeding. Any bleeding or tenderness which may occur initially may soon fade as the gum tissues become firm and healthy. (See your Dentist if this persists.)

**Sulcabrush**<sup>®</sup>

Replacement Tips

**POCKET™ Sulcabrush**<sup>®</sup>  
with 1 extra  
Replacement Tip

POCKET Pack  
**Sulcabrush**<sup>®</sup> with  
**POCKET™ Sulcabrush**<sup>®</sup>

**Travel Pouch**